

XENLING



FIELD GUIDE / XL-013
LIVE CULTURE

SETUP + CARE

UNFOLD FULLY

Keep this guide.

SUBSTRATE ALREADY INSIDE?



Alien containment
unit

Alien substrate already in your
alien containment unit?
Skip Step 0. Turn over to Step 1.
Powder kit? Start below.

0a MIX / POWDER KIT ONLY



0.5 g



50 ml

Mix all 0.5 g alien substrate
powder with 50 ml water
in a small pan.
50 ml makes a batch, not a dose.

0b BOIL + COOL



Bring to a boil; stir until dissolved.
Remove from heat. Adult prep.
For pipetting: cool until warm,
still liquid. Never use it boiling.

0c FILL THE EMPTY UNIT



RATED?



Heat-rated 0.3 ml pipette only:
add a shallow layer while warm.
Rating unknown? Let substrate
set + cool; use clean tweezers
to add a small piece instead.

0d LET IT SET



Alien substrate must set and cool
to room temperature before it
touches your Xenling. If the unit
already holds a live sample,
add only fully cooled pieces.

1 SETTLE IN



Keep upright at room temperature, away from direct sunlight.
Your Xenling needs air:
do not seal the unit airtight.

2 FEED YOUR XENLING



Place a small piece of alien food beside your Xenling.
Moisten food with clean water,
one drop at a time.

3 KEEP THE SUBSTRATE DAMP



If alien substrate starts drying,
add clean, room-temperature
water one drop at a time.
Damp, never flooded.

4 WATCH IT EXPLORE



Check daily for branching growth.
Offer a small fresh piece of alien
food as it grows. Feed sparingly.
Remove old, uneaten food.

5 CLEAN + RESET



Fuzz or bad smell?
Remove old alien food.
If it persists, save a small,
healthy alien sample.

Clean + rinse the unit.
Add fresh, cooled alien substrate.
Return your saved sample.
**Contact seller for alien
substrate powder.**

HANDLE WITH CARE



**DO NOT
RELEASE**

Wash hands before + after care.
Do not eat. Do not release.
Pipette: 0.3 ml capacity,
not a fixed watering dose.